
药剂师能看到的情况越多，那么他改进某一流程就越方便
The more the chemist can see, the easier it is for him to influence or change a process

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Glass plants in kilolabs

Glass is the ideal material for plants producing kilograms of product required for trials

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A kilolab is usually employed at some stage between a research laboratory and a pilot plant or production facility. As the name implies, its main purpose is to produce kilo amounts of products. These are often products that have come out of a research laboratory and have not been tested in a stirred vessel; so many parameters are still unknown. Kilolabs also serve the function of developing processes that can be easily duplicated in the pilot plant or in a production facility.

即使两个部分不是完全的对准，每个接头也都能由四颗带有弹簧的螺钉很好的固定
Assembling the join using four screws with springs gives a perfect seal, even if the two pieces are misaligned

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The chemist who gets a receipt from the research laboratory has to quickly produce considerable amounts of product. He or she needs a flexible unit, but it should not be fully automated, or if it is it should have a local override.

Such a unit consists usually of a reactor with a working volume of 15 to 250 litres. This is an effective way of stirring the content, usually with a retreat curve impeller and a range of about 0 to 350 revolutions per minute. A heating and cooling unit is also required, whose range depends on the processes.

Solids are directly added through an opening on the reactor. Liquids are either added through the feed vessels by opening valves or by dosing pumps, which can accurately change the flow rates.

In most cases, equipment is needed to distil and separate phases. Distilling is usually done with one or several main condensers in parallel and a vent condenser to protect the vacuum pump from any chemicals. The phases can be separated in the reactor through the bottom outlet valve, or with the distillate through a phase separator. At this point sometimes columns are needed.

The distillate is received in a secondary reactor or a simple glass vessel to collect solvents. Condensers will need cooling, using tap water, an in house glycol system or a separate chiller.

A nitrogen source is needed to overlay the system with the inert gas, and a vacuum source needs to be hooked up for the distillation process. The exhaust gases during the process need to be scrubbed before they can be released into the environment. This can be done with a dedicated scrubber in glass. Filters are sometimes needed for further processes, or other equipment such as centrifuges.

A kilolab can also be used to refine, change or re-evaluate a process. It can also be assigned to a production facility to quickly solve problems showing up during production runs.

Why use glass?

The most obvious critical parameter is chemical resistance. Glass is resistant to all acidic reactions and some of the alkaline reactions. Borosilicate glass is an obvious choice, since units have to be used for a wide range of processes that have already been conducted in a research lab where borosilicate glass has been used. As well as having the advantage of transparency, glass is less costly than any metal.

The more the chemist can see, the easier it is to influence or change a process. In many cases the chemist likes to be able to intervene directly; that is, manually. He or she also needs to record and evaluate what is happening. This means parameters including temperature, pressure, agitator speed and pH.

Because of these advantages, glass is the material of choice for kilolabs.

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
Safety issues

Using glass is not a problem in the laboratory, but with bigger equipment safety departments are hesitant to use glass as a material for process equipment. In Europe there are strict regulations and calculations available to obey before using glass for this purpose. This is why these plants are limited in pressure from full vacuum to a maximum of 0.5 bar, and all units are further protected with rupture discs. Sometimes the units are installed into hoods to further protect the operator, but also to limit the environment which needs to be explosion proof.

Heating and cooling

Another critical parameter in the choice of material is its ability to be heated and cooled. Glass is a good insulator, a property that improves with thickness. That's why the best reactors are glass-lined steel vessels in the bottom section, with glass above the working volume – again for visual reasons. This allows a temperature range of -60°C to +200°C, or as low as - 80°C if

使用玻璃对实验室来说没有问题，不过安全部门在使用玻璃制造一些更大型的仪器时还是比较犹豫
Using glass is not a problem in the laboratory, but with bigger equipment safety departments are hesitant to use glass as a material for process equipment

玻璃可以耐受所有的酸性反应和部分碱性反应
Glass is resistant to all acidic reactions and some of the alkaline reactions

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Protecting a kilolab against explosions

Chemists at kilolabs deal with large amounts of chemicals and solvents. Even if they are contained, everything in the same room or hood needs to be explosion-proof. Typical parameters are:

Temperature: The reactor temperature needs to be controlled. Depending on the heating cooling unit chosen the control might come with such units. The local indication is easily solved with an explosion proof indicating transmitter, whereas the change of set point needs an explosion proof panel with a set point indication and a knob for adjustment. Other temperatures to at least monitor are vapour and condensate temperatures and if there is a column installed each stage needs another temperature measurement.

Pressure: Since the reactions are done at atmospheric pressure and during a distillation as much of a vacuum is drawn as possible sometimes a simple indication with a gauge is sufficient. For a recording of the pressure a sensor is necessary with a local transmitter and indication.

Agitator speed: The changing of the rpm can be done with a stepless mechanical transmission or a variable frequency drive. Control is not so much an issue since rpm is not a critical measurement but should be known. Also it is rarely changed during a process.

Other parameters are Flow rates and pH

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
special glass lining is used). For reactions with temperatures down to -125°C, a complete Hastelloy vessel without a glass section is needed, with a specially designed jacket for liquid nitrogen. In all these cases the glass upper section can remain glass, since all separations will occur at somewhat elevated temperatures. The glass upper section is designed for a range of -60°C to +200°C. This limitation is not given by the glass but is dictated by the use of PTFE for the gaskets.

Vacuum tightness

The vacuum tightness of the whole plant is another critical parameter – and not just to minimise the size of vacuum pump required. As introducing oxygen is usually detrimental to the process, it is necessary to be able to create an inert environment. This can be achieved by overlaying nitrogen and minimising the leak rate at each joint or seal.

Büchi discovered a unique way of doing this by designing each joint with a ball and socket arrangement. The socket is a perfectly flame polished sphere, whereas

除了拥有透明的特点，玻璃与金属材料相比还有价格便宜的优势
As well as having the advantage of transparency, glass is less costly than any metal

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the ball has a machined PTFE gasket. Assembling the join using four screws with springs gives an almost perfect seal, even if the two pieces are misaligned. The four springs will compensate for any misalignment of several degrees, without compromising vacuum tightness and without requiring any additional flexible elements. Because PTFE flows, especially at elevated temperatures, the springs will also compensate for any movement. All this together will assure a minimum leakage rate and reduce the amount of oxygen introduced into the reactor environment. The leak rate of a medium-sized plant using this system is well below 5 mbar per hour. ■

最好的反应器是下部使用玻璃线纹的钢管，工作体积上方配以玻璃
The best reactors are glass-lined steel vessels in the bottom section, with glass above the working volume

工艺设备 Process Equipment

3D engineering

Large glass overheads

distillation

API's

synthesis

rectification

crystallization

extraction

中试反应设备 Pilot Plants and kilo Lab

minipilot
5 to 10 liters

chemreactors
15 to 250 liters

Nutsche-Filters
10 to 50 liters

Scrubbers
10 to 50 liters

**büchiflex –
the unique flexible
glass connection**

safety & performance

安全的压力反应 safe pressure reactions

synthesis

**Change the
vessel – not the
entire system!**

polymerization
catalyst testing

hydrogenation
crystallization

**tinyclave
miniclave drive**
10 – 25/50 – 300 ml

Vessel type 1

Borosilicate glass
max. 12 bar

Vessel type 3

stainless steel,
Hastelloy, Tantalum,
Titanium, Zirkonium

picoclave
50 – 300 ml

polyclave
0,25 – 5 l

kiloclave
0,5 – 20 l

BuchiGlasUster China Co., Ltd.
比棋优仪器设备(上海)有限公司
上海市徐汇区肇家浜路1065号1406室
Room 1406, No. 1065, Zhaojiafang Road.
Xuhui District, Shanghai 200030, China

TEL: (021) 6468-6338, FAX: (021) 6468-8317
手机: 13900699339, 15921773386
邮箱: mail@bguchina.com or william.wjh@gmail.com
网址: www.bguchina.com

 **büchiglasuster**
switzerland
Pilot Plants and Reactor Systems
www.buchiglas.ch