

# **JULABO PRESTO® W91tt**

Temperature stability with a 100 l reactor at -50 °C

## **Objective**

This case study tests the temperature stability of a JULABO PRESTO<sup>®</sup> W91tt with a 100 liters glass reactor. The W91tt is connected to the reactor via two 2.0 m metal tubings. The W91tt was set to a set point of -50 °C.

JULABO PRESTO® W91tt

# **Test Conditions**

JULABO unit Cooling power

Heating capacity Band limit Flow pressure Bath fluid Reactor

Control

+20 °C | 11.0 kW 0 °C | 10.0 kW -20 °C | 9.5 kW 36 kW ohne 0.45 bar JULABO Thermal HL80 100 liters glass reactor (Büchiglas) filled with 100 liters Thermal HL80 External (ICC)

#### **Test Results**

See chart on back page: The W91tt cools down the reactor to -50 °C. After reaching the temperature of -50 °C, the temperature within the reactor fluctuated by  $\pm 0.01$  K max.



Hz

#### Environment

Room temperature	+20 °C
Humidity	45 %
Voltage	3 x 400 V / 50



**Tip** You can also use the robust Pt100 with PTFE coating.


More tips on back page >>



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



#### www.julabo.de



#### Tip

Make use of the option to regulate the pump pressure. You can define the desired pressure in the PRESTO® settings.



**Tip** The Ethernet interface permits full access to all operational functions of the PRESTO®.



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



## www.julabo.de