

JULABO PRESTO® W92tt

Temperature stability with a 100 l reactor at -50 °C

Objective

This case study tests the temperature stability of a JULABO PRESTO[®] W92tt with a 100 liters glass reactor. The W92tt is connected to the reactor via two 2.0 m metal tubings. The W92tt was set to a set point of -50 °C.

JULABO PRESTO® W92tt

Test Conditions

JULABO unit Cooling power

Heating capacity Band limit Flow pressure Bath fluid Reactor

Control

+20 °C | 19.0 kW 0 °C | 15.5 kW -20 °C | 9.5 kW 36 kW ohne 0.45 bar JULABO Thermal HL80 100 liters glass reactor (Büchiglas) filled with 100 liters Thermal HL80 External (ICC)

Test Results

See chart on back page: The W92tt cools down the reactor to -50 °C. After reaching the temperature of -50 °C, the temperature within the reactor fluctuated by ± 0.02 K max.



Environment

Room temperature	+20 °C
Humidity	45 %
Voltage	3 x 400 V / 50 Hz



Tip You can also use the robust Pt100 with PTFE coating.


More tips on back page >>



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



www.julabo.de



Tip

Make use of the option to regulate the pump pressure. You can define the desired pressure in the PRESTO® settings.



Tip The Ethernet interface permits full access to all operational functions of the PRESTO[®].



JULABO GmbH Eisenbahnstraße 45 77960 Seelbach / Germany Tel. +49 (0) 7823 51-0



www.julabo.de