

PRESTO[®] W92tt Cooling a 100 liters reactor

from +20 °C to -80 °C

Objective

0

This case study tests the cooling power of PRESTO[®] W92tt with a 100 liters glass reactor. The PRESTO[®] W92tt is connected to the reactor via two 3 m metal tubings. The PRESTO[®] W92tt is programmed to cool down from +20 °C to -80 °C.

Environment

Room temperature	+20 °C
Humidity	45 %
Voltage	400 V / 50 Hz

Test Conditions

JULABO unit Cooling power

Heating capacity Band limit Flow pressure Bath fluid Reactor

Jacket volume Control +20 °C 19 kW 0 °C 15.5 kW -20 °C 9.5 kW 36 kW with 0.5 bar Thermal HL80 100 I glass reactor (Büchiglas) filled with 70 I Ethanol 30 I External (ICC)

PRESTO® W92tt








Test Results

0

The PRESTO® W92tt cooling process from +20 °C to -80 °C in 3 h 45 min without overshoot.



Measured with EasyTEMP Professional

Tip

5

Take advantage of our wide range of accessories. The M+R adapter enables you to display and record an additional temperature.



Tip

You can also use the robust Pt100 with PTFE coating.



h