

PRESTO® W92tt

Cooling a 100 liters reactor from +20 °C to -20 °C

Objective

This case study tests the cooling power of PRESTO® W92tt with a 100 liters glass reactor. The PRESTO® W92tt is connected to the reactor via two 3 m metal tubings. The PRESTO® W92tt is programmed to cool down from +20 °C to -20 °C.

Environment

Room temperature +20 °C 45 % Humidity

400 V / 50 Hz Voltage



Test Conditions

JULABO unit PRESTO® W92tt Cooling power +20 °C 19 kW 0 °C 15.5 kW

-20 °C 9.5 kW

Heating capacity 36 kW **Band limit** with Flow pressure 0.5 bar Bath fluid Thermal HL80

Reactor 100 l glass reactor (Büchiglas)

filled with 70 I Ethanol

Jacket volume 30 I

Control External (ICC)








Test Results

The PRESTO® W92tt cooling process from +20 °C to -20 °C in 1 h 30 min without overshoot.



Setpoint

Temperature in reactor's interior

Temperature in reactor's jacket

Tip

Take advantage of our wide range of accessories. The M+R adapter enables you to display and record an additional temperature.



Tip

Use the free of charge EasyTEMP software to control the units with the PC and to show the temperature curves graphically.

