

PRESTO® W92tt

Heating a 100 liters reactor from 0 °C to +20 °C

Objective

This case study tests the heating power of PRESTO® W92tt with a 100 liters glass reactor. The PRESTO® W92tt is connected to the reactor via two 3 m metal tubings. The PRESTO® W92tt is programmed to heat up from 0 °C to +20 °C.



Room temperature +20 °C 45% Humidity

400 V / 50 Hz Voltage



Test Conditions

JULABO unit PRESTO® W92tt Cooling power +20 °C 19 kW 0 °C 15.5 kW

-20 °C 9.5 kW

36 kW

Heating capacity **Band limit** with Flow pressure 0.5 bar Bath fluid Thermal HL80

Reactor 100 liters glass reactor (Büchiglas)

filled with 70 I Ethanol

Jacket volume 30 I

Control External (ICC)








Test Results

The PRESTO® W92tt heating process from 0 °C to +20°C in 35 min without overshoot.



Setpoint

Temperature in reactor's interior

Temperature in reactor's jacket

