

## CAT Stirrer Motors



**CAT R100 Series** - high torque stirrer motors in four different executions for stirring volumes up to 200 liters of water.

**CAT R20** - lightweight, maintenance-free stirrer motor for volumes up to 30 liters of water.

Article No.	Type	Engine output	Revolutions Rpm	Torque Ncm	Chuck mm	Weight kg
23.658.10	R20	60 W	50 - 2000	20	0.5 - 8	1.3
23.658.20	R100 SD	300 W	40 - 2000	100	1 - 10,5	5.0
23.658.30	R100 C	300 W	40 - 2000	150	1 - 10,5	5.0
23.658.40	R100 CL	300 W	30 - 1100	250	1 - 10,5	5.0
23.658.50	R100 CT	300 W	10 - 500	600	1 - 10,5	5.5



*CAT stirrer motors stand out with an excellent price-performance ratio. Complete technical specifications you can find here: <http://www.cat-ing.de>*

*We also can offer stirrer motors from other manufacturers.*